



Consumable pricing EasyOne Pro

Validity: one Year (until April 14th. 2012)
 Effective: April 14th. 2011
 Conditions: net ex works, no discount applicable

Consumable:




Picture	Product name	Order #	Price (Contact your Dealer)	Extras
	Spirette box of 50	2050-1		
	DLCO Bariette box of 50	3050-1		

Note:










Spirette: 1 for each person to be tested, single patient use is recommended.
 For prices refer to our price list on our homepage (dealer access on our web)

Barriette: 1 for each person to be tested, single patient use is recommended.
 This is the part that covers the valve and protects it and the EasyOne Pro from contamination by patients: barriette is a filter.

Recommended change parts (once a year):

Picture	Product name	Order #	NET Sales Price (CHF)	Extras
	Filter pack	3000-50.1sp		
	Gas supply tube incl. filter, sample draw tube, conditioning unit	3000-50.8sp		
	One-way valve	3000-50.9sp		

Information for selected spare parts:

Picture	Product name	Order #	NET Sales Price (CHF)	Extras
	Filter pack	3000-50.1sp		
	Motor block	3000-50.2sp		
	DLCO Valve	3000-50.3sp		
	Pressure tube (FDA approved)	3000-50.4sp		
	Quick Connector for Gas Regulator	3000-50.4.1sp		
	Flowsensor	3000-50.5sp		
	Stylus touch	3000-50.6sp		
	Temperature sensor	3000-50.7sp		
	Gas supply tube incl. filter, sample draw tube, conditioning unit	3000-50.8sp		
	One-way valve	3000-50.9sp		
	Overflow Valve	3000-50.10sp		